

Instrumentation Designed and Manufactured by Labfacility

We offer a large range of temperature instrumentation including our flagship range of bench top units, designed and manufactured by Labfacility. These are sold Worldwide directly and through our distributors and are suitable for a large number of applications.

L60+ Thermocouple & Fine Wire Welder

The Thermocouple Welder is a compact, simple-to-use instrument designed for thermocouple and fine wire welding. Use of the Thermocouple Welder does not require special skills and most operatives will be capable of producing quality work with little practice. The instrument is supplied with a full range of user accessories including a footswitch. [More details.](#)



Tempmaster PRO Precision Thermometer

The Tempmaster PRO from Labfacility provides very high precision temperature measurement with Pt100 sensors. Semi Standard Probe and calibration are also available with this product. [More details.](#)

Labcal PRO Precision Thermometer

The Labcal PRO from Labfacility is fully characterised for Pt100 sensors and all major thermocouples. As with the Tempmaster, Semi Standard Probe and calibration are available with this product. [More details.](#)



L200 - 8 Channel Temperature Monitor Pt or TC inputs

The Labfacility L200-PT for Pt100 & L200-TC Thermocouple Thermometers (Temperature Loggers) can be used in conjunction with a PC to provide accurate, versatile 8 channel thermocouple temperature measurement, scanning and logging of measured values. Bench Top or Panel Mount options are available for the L200 range. [More details.](#)

Product videos are available to view on our YouTube channel.