


4-Kern-PFA-Isolier, verzinnbeschichteter Kupferschirm, PRT-Sensorkabel / -kabel

10M




4 CORE
XE-2734-010

25M



4 CORE
XE-2734-025

10M



4 CORE
XE-2734-010

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

4-Kern-PFA-Isolier, verzinnbeschichteter Kupferschirm, PRT-Sensorkabel / -kabel

Hohe Temperatur bis zu 260C, extrudierte PFA-Teflonkonstruktion. Hufig fr Pt100-Thermometer verwendet, auch geeignet fr andere Typen wie Pt130, Pt500 und Pt1000

Technische Daten**Specifications**

General Description	PFA insulated 7/0.2mm tin plated copper conductors, twisted cores, tin plated copper braided screen with PFA insulated outer jacket
Cable Length	10m, 25m, 50m, 100m
Cable Type	PFA Teflon® construction
Core / Strands	7/0.2mm
Insulation Material	PFA Insulated
Max. Temperature	+260°C
Min. Temperature	-75°C
Number of Wires	4
Standards Met	Cores colour coded Red & White in accordance with IEC-751, black outer jacket