



## 4-Kern-PTFE-Isolierter, versilberter Kupferschirm, PRT-Sensorkabel / -leitung



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**4-Kern-PTFE-Isolierter, versilberter Kupferschirm, PRT-Sensorkabel / -leitung**

Kann verwendet werden, um 2-, 3- oder 4-adrige RTD-Konfigurationen zu erweitern. Hufig für Pt100-Thermometer verwendet, auch geeignet für andere Typen wie Pt130, Pt500 und Pt1000

**Technische Daten****Specifications**

<b>General Description</b>	4 core PTFE insulated 7/0.2mm silver plated copper conductors, twisted cores, kapton tape, silver plated copper braided screen with PTFE insulated outer jacket
<b>Cable Length</b>	10m, 25m, 50m, 100m
<b>Cable Type</b>	wrapped PTFE tape construction
<b>Core / Strands</b>	7/0.2mm
<b>Insulation Material</b>	Overall black PTFE insulation
<b>Max. Temperature</b>	+260°C
<b>Min. Temperature</b>	-75°C
<b>Number of Wires</b>	4
<b>Standards Met</b>	Cores colour coded 2 x Red / 2 x White to IEC-751 + black outer jacket