



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PTFE-isoliertes, doppelt verdrehtes IEC-Thermoelement mit freiliegender Sperrschicht und Miniaturstecker - Typ T - 300 mm Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Thermoelement Typ K, J oder T. Schnell ansprechende, geschweifte, freiliegende Verbindungsstelle, erhältlich in der Kabelspezifikation der Klasse 1 gem BS EN 60584-3:2008 mit einem Miniatursteckeranschluss.

Ideal für Test- und Entwicklungsanwendungen.

Leistungsbeschreibung

Specifications

Product Code	XE-0452-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	T IEC
Cable Length	300mm
Cable Type	Teflon® PTFE insulated twin-twisted lead
Core / Strands	1/0.2mm
Termination Type	Miniature flat pin plug colour coded in accordance IEC 584
Max. Temperature	260°C
Min. Temperature	-75°C
Tolerance	Class 1
Standards Met	IEC-584-3
Connector Type	Plug
Connector Size	Miniature
Series	Exposed Junctions