




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

ANSI PFA isoliertes Thermoelement mit freiliegender Verbindungsstelle und blanken Schwänzen - Typ K - 1/0,3 mm Stranggare - 2 m Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-3544-001
General Description	Fast response, welded exposed junction
Sensor Type	Exposed Junction
Thermocouple Type	k
Cable Length	2M
Cable Type	Teflon® PFA insulated twin-twisted lead
Core / Strands	1/0.3mm
Termination Type	Bare Tails
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	ANSI MC96.1
Series	Exposed Junctions