




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

ANSI Glasfaserisoliertes Thermoelement mit Sichtschutz und blanken Schwnzen - Typ K - 2M Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-3546-001
General Description	Fast response, welded exposed junction
Sensor Type	Good temperature resistance (not suitable in fluids)
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	2M
Cable Type	glassfibre insulated flat pair cable
Core / Strands	1/0.3mm
Cable Termination	Bare tails
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	ANSI MC96.1