




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

Selbstklebender Patch PT100 Sensor - Typ PRT - 1M Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-3662-001
General Description	Pt1000 element encapsulated in a silicone rubber patch, Self-adhesive, aluminium foil backing
Sensor Type	Pt1000 (1000 Ω @ 0°C), thin film Class B to IEC751
Element	Pt1000 sensor to IEC751 Class B
Patch dimensions	Length 30mm x Width 15mm x Height 4mm
Cable Length	1M
Cable Type	Teflon ® insulated twin twisted
Core / Strands	7/0.2mm
Cable Termination	Bare tails
Max. Temperature	150°C
Min. Temperature	-50°C
Series	Patch Pt1000