


## **Lufttemperatur / PT100 Sensor fr den Innenbereich**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

<b>Product Code</b>	XE-3670-001
<b>General Description</b>	PT100 wall mounting indoor temperature sensor for environmental temperature measurement
<b>Sensor Type</b>	For indoor use only
<b>Housing</b>	Moulded case L85 x W85 x H30mm
<b>Max. Temperature</b>	+70°C
<b>Min. Temperature</b>	-50°C
<b>Accuracy</b>	Class B to IEC 751
<b>Number of Wires</b>	Can be connected to any 2, 3 or 4 wire instrument
<b>Features #1</b>	Vented construction
<b>Features #2</b>	Surface or standard electrical conduit mounting
<b>Features #3</b>	Cable entry slots underside of housing
<b>Series</b>	Air Indoor Sensor