



Handgehaltenes federbelastetes Thermoelement mit Kupferscheibenspitze Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Metern gewickeltem Verlängerungskabel, abgeschlossen in einem Miniatur-Thermoelement-Stecker.

Specifications

Product Code	XE-3816-001
General Description	Surface temperature spring loaded thermocouple with copper disc tip
Thermocouple Type	K
Calibration	IEC
Cable Length	2m Fully Extended
Cable Type	Retractable Curly Thermocouple Lead - Polyurethane outer
Cable Termination	Miniature Plug
Max. Temperature	+600°C
Dimensions	Probe 4.7mmO x 63mm long stem, O8 x 17mm tip with O4.5mm copper disc
Tip	Copper disc tip
Series	Handheld Probe