


## Typ K IEC rechtwinklig Kupferscheibenspitze



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Metern gewickeltem Verlängerungskabel, abgeschlossen in einem Miniatur-Thermoelement-Stecker.

## Specifications

<b>Product Code</b>	XE-3817-001
<b>General Description</b>	Surface temperature spring loaded thermocouple with right angled probe and copper disc tip
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+600°C
<b>Dimensions</b>	Probe O8mm(17mmL) Tip x 60mm Hot Leg and O6mm x 180mm Cold Leg with O4.5mm copper disc.
<b>Series</b>	Handheld Probe