

Oberflächentemperatur-Sonde Typ K IEC 90 mm lang mit 10 mm Spitze



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Metern gewickeltem Verlängerungskabel, abgeschlossen in einem Miniatur-Thermoelement-Stecker.

Specifications

Product Code	XE-3818-001
General Description	Surface temperature fast response thermocouple with ceramic tip and coiled element
Thermocouple Type	K
Calibration	IEC
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	+900°C
Dimensions	Probe O6 x 90mm long stem with O12 x 10mm tip. Coiled element diameter 7mm
Series	Handheld Probe