


## Typ K IEC rechtwinklige Oberfläche



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Edelstahlmantel mit geformtem Griff und 2 Metern gewickeltem Verlängerungskabel, abgeschlossen in einem Miniatur-Thermoelement-Stecker.

## Specifications

<b>Product Code</b>	XE-3819-001
<b>General Description</b>	Surface temperature fast response thermocouple
<b>Sensor Type</b>	Right angled Probe
<b>Thermocouple Type</b>	K
<b>Calibration</b>	IEC
<b>Cable Length</b>	2m
<b>Cable Termination</b>	Fitted Plug
<b>Max. Temperature</b>	+900°C
<b>Dimensions</b>	Probe O12mm Tip x 40mm Hot Leg and O6mm x 180mm Cold Leg. Coiled element diameter 7mm
<b>Tip</b>	Ceramic Tip, Coiled Element
<b>Series</b>	Handheld Probe