



100 Ohm Flachschiichtdetektoren - Klasse 1/3 DIN - 2mm x 5mm



Labfacility are currently multi-sourcing flat film detectors in order to ensure we can maintain a supply service to our customers without affecting the quality or functionality for the foreseeable future.

All flat film detectors supplied, although differing in appearance and packaging in some cases, will meet the specifications as published on our data sheets.



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Labfacility beschafft derzeit Flachsichtdetektoren mehrfach, um sicherzustellen, dass wir unseren Kunden auf absehbare Zeit einen Lieferservice bieten können, ohne die Qualität oder Funktionalität zu beeinträchtigen.

Alle gelieferten Flachfilmdetektoren unterscheiden sich zwar teilweise in Aussehen und Verpackung, erfüllen aber die Spezifikationen, die auf unseren Datenblättern veröffentlicht sind.

Leistungsbeschreibung

Specifications

Product Code	XE-4460-001
Sensor Type	Pt100 (100 Ohms @ 0°C)
Element	Thin film, 10mm tails
Max. Temperature	500°C
Min. Temperature	-50°C
Tolerance	1/3DIN in accordance with IEC 60751
Dimensions	2mm x 5mm
Fundamental interval	38.5 ohms (nominal)
Ice point resistance	100 ohms
Self-heating	<0.5°C/mW
Stability	±0.05%
Thermal response	<0.4s @ t0.9 (in water at 0.4m/sec)
Important Information	Resists vibration & shock
Series	RTD Detectors