




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

100 Ohm drahtgewickelte Detektoren - Klasse A - 2,8 mm x 15 mm



100 Ohm drahtgewickelte Detektoren - Klasse A -
2,8 mm x 15 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bitte beachten Sie: Die Lieferzeiten der Klasse 1/10 DIN variieren - bitte erkundigen Sie sich oder Sie werden mit einer Vorlaufzeit kontaktiert, nachdem Sie Ihre Bestellung aufgegeben haben.

Leistungsbeschreibung

Specifications

Product Code	XE-4746-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	Class A
Dimensions	2.8mm Width x 15mm Length
Series	RTD Detectors