




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

100-Ohm-Draht-Detektoren Klasse 1/10 DIN 1,5 mm x 8 mm



100 OHM 1/10 DIN 1,5 x 8 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Bitte beachten Sie: Die DIN-Lieferzeiten der Klasse 1/10 können variieren bitte fragen Sie nach, damit wir Sie vor der Bestellung informieren können.

Technische Daten

Specifications

Product Code	XE-4788-001
General Description	Wirewound construction for optimum performance and stability
Sensor Type	Pt100 Single Element
Max. Temperature	650°C
Min. Temperature	-200°C
Tolerance	1/10DIN relative to IEC 60751 over the range 0°C to 150°C
Dimensions	1.5mm Width x 8mm Length