

## GlimmerdÃfÂ¼senheizungen



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Diese Produkte sind verfügbar, solange der Vorrat reicht, mit begrenzten Mengen für jede Version, Sie werden benachrichtigt, wenn Sie in den Warenkorb gelegt werden, wenn wir nicht genügend Lager haben - [bitte kontaktieren Sie uns](#), wenn Sie spezifische Anforderungen haben.

## Specifications

|                              |   |
|------------------------------|---|
| <b>General Description</b>   | Best quality Mica for high dielectric strength and maximum heat transfer  |
| <b>Element</b>               | High performance precision wound element for superior heat distribution   |
| <b>Cable Length</b>          | 1 metre long stainless steel braided cable  |
| <b>Max. Temperature</b>      | Recommended for temperatures up to 300°C  |
| <b>Output</b>                | 240V  |
| <b>Features #1</b>           | Full width low expansion straps   |
| <b>Features #2</b>           | Specially treated finish to resist oxidation and enhance thermal conductivity   |
| <b>Important Information</b> | Wide range of applications including plastic process machinery, injection / blow moulding, rubber moulding, test equipment, food process industry etc |