


Magnetband-Thermoelement mit Miniaturstecker Typ K IEC - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Ein flexibles Thermoelement mit Magnetstreifen, das an den meisten ferromagnetischen Oberflächen befestigt werden kann, die große Fläche des Magneten sorgt für eine gute Haftung und ermöglicht eine einfache Installation. Nützlich in einer Vielzahl von Anwendungen zur Temperaturüberwachung, insbesondere in Test und Entwicklung.

Leistungsbeschreibung

Specifications

Product Code	XE-3455-001
General Description	Ideal for ferrous metal surface measurement
Thermocouple Type	K
Thermocouple Junction	Grounded measuring junction
Magnet	50 x 25 x 1mm magnetic strip
Cable Length	2M
Cable Type	PTFE insulated wires, twin twisted
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	100°C
Min. Temperature	-50°C
Standards Met	IEC
Series	Magnet Thermocouple