



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC Rohroberfläche Jubilee Clip Style Thermoelement, Glasfaser-Edelstahlgeflecht - Typ K - passend fr 30-45 mm Rohrdurchmesser



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-3916-001
General Description	Stainless steel pipe clamp sensor, jubilee clip style, with adjustable diameter
Thermocouple Type	K
Thermocouple Junction	Grounded tip for fast response
Cable Length	2M
Cable Type	Glassfibre stainless steel braided extension cable
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	350°C
Standards Met	IEC
To Fit Pipe Diameters	30-45mm