



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Glasfaserisoliertes flaches Paar IEC-Thermoelement mit freiliegender Sperrschicht - Typ K mit Bare-Tails-Anschluss - 1/0,508 mm x 5 m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5534-001
General Description	Good temperature resistance (not suitable in fluids)
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	5M
Cable Type	Glassfibre
Core / Strands	1/0.508mm
Termination Type	Bare Tails
Max. Temperature	+350°C (short periods up to 400°C)
Min. Temperature	-60°C
Tolerance	Class 2
Standards Met	IEC-584-3
Series	Exposed Junctions