




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Platin-Widerstandsdetektoren mit verlängerten Leitungen, drahtgewickelte Keramik - 2 Drhte - 300 mm Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Optimale Leistung, geeignet fr Luft- und Tauchanwendungen (Geeignet nur in trockenen, nicht kontaminierenden Atmosphren, wenn sie ungeschutzt sind. Geeigneter Schutz / Ummantelung erforderlich fr das Eintauchen in Flligkeiten, korrosive Substanzen, Prozesse usw.).

Leistungsbeschreibung

Specifications

Product Code	XE-5595-001
General Description	Platinum Resistance Detectors with extended leads, Wire-Wound Ceramic
Sensor Type	Platinum Resistance, Pt100 Ohm
Sheath Material	Platinum coil wire-wound construction sealed inside a high purity alumina ceramic body
Cable Length	300mm
Cable Type	Glass-fibre insulated nickel
Cable Termination	Bare tails
Max. Temperature	350°C
Min. Temperature	-50°C
Standards Met	European / British Standard Conforms to IEC 60751:2008 / BS EN 60751:2008
Accuracy	Class B
Extended lead wires	2 Wire
Detector Dimensions	2.8mm x 15mm
Fundamental interval	38.5 ohms nominal (100 ohms)
Ice point resistance	100 Ohms
Series	Wire Wound Detector with Extended Leads