



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

PFA-Kabel Tidy IEC Thermoelement mit freiliegender Sperrschicht und montiertem Ministecker - Typ K - PFA Flat Pair - 1M Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Kabelgebundenes Thermoelement mit eingebautem Miniatur-IEC-Stecker, PFA Twin Twisted oder Flat Pair.

Ein hochempfindliches Präzisions-Thermoelement mit feinem Draht und offener Verbindung mit einem passenden Miniaturstecker.

Leistungsbeschreibung

Specifications

Product Code	XE-5599-001
General Description	A complete, precision fine wire, high sensitivity thermocouple with a fitted miniature connector and cable tidy.
Sensor Type	Exposed Junction
Thermocouple Type	K IEC
Cable Length	1M
Cable Type	PFA Flat Pair insulated wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	Rated 220°C continuous
Tolerance	Class 1
Standards Met	IEC-584
Connector Type	Plug
Connector Size	Miniature
Series	Cable Tidy Exposed Junction