


IEC Twist Melt Bolt Thermoelement - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5532-001
General Description	Rotational blade style tip allowing correct positioning for reduced resistance in the melt flow
Sensor Type	Widely used in plastics extruders and injection moulders
Thermocouple Type	K
Probe Length	152mm stem length
Sheath Material	Stainless steel bolt with 1/2"UNF-20 2A thread
Termination Type	Standard plug termination (round pin)
Max. Temperature	400°C
Tolerance	Tolerance & colour code in accordance with IEC-584 specification
Standards Met	IEC
Thread	1/2"UNF-20 2A thread
Tip	20mm long blade tip immersion length