

IEC Velcro Thermoelement, PFA mit Miniaturstecker - Typ T



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5611-001
General Description	Velcro thermocouples can be used as a quick and simple method of measuring the temperature around the surface of a pipe.
Thermocouple Type	T
Cable Length	1.5M (0.5M Inside Velcro)
Cable Type	PFA
Core / Strands	7/0.2mm
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Max. Temperature	150°C
Standards Met	IEC
Length	Velcro: 500mm