



LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

IEC-Feindraht-Scheiben-Thermoelement, PTFE mit Miniaturstecker - Typ K - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5621-001
General Description	Exposed junction thermocouple with Copper Disc brazed at tip
Sensor Type	Copper disc
Thermocouple Type	Type K
Thermocouple Junction	Grounded measuring junction
Cable Length	2M
Cable Type	1/0.2mm PTFE insulated twisted wire
Termination Type	Miniature flat pin plug, colour coded in accordance IEC 584
Standards Met	IEC
Series	Copper Disk Thermocouple