


ANSI Hochtemperatur-Glasfaser-Thermoelement mit zwei verdrehten, freiliegenden, geschweiften Spitzen Typ K mit Standardsteckeranschluss - 20 m Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5632-001
General Description	Fast response, welded exposed junction. Good temperature resistance (not suitable in fluids)
Sensor Type	Exposed Junction
Thermocouple Type	K
Cable Length	20M
Cable Type	high temperature glassfibre insulated twisted cable
Core / Strands	1/0.711mm conductors (21awg)
Cable Termination	Standard Plug
Max. Temperature	700°C
Min. Temperature	0°C
Tolerance	Tolerance to ANSI MC96.1 Class 2
Standards Met	Colour code to ANSI MC96.1
Connector Type	Plug
Connector Size	Standard
Colour	Yellow
Series	Exposed Junctions