




LABFACILITY
TEMPERATURE & PROCESS TECHNOLOGY

PFA-isoliertes ANSI-Thermoelement mit freiliegender Sperrschicht und angespritztem Stecker - Typ K - 2M



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

Bitte beachten Sie, dass Sie bezüglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.

Leistungsbeschreibung**Specifications**

Product Code	XE-0026-001
General Description	Good mechanical strength and flexibility, resistant to oils, acids and other adverse fluids, gasses, steam etc
Thermocouple Type	K
Thermocouple Junction	Exposed
Cable Length	2M
Cable Type	PFA Teflon® twin twisted pair
Core / Strands	1/0.3mm
Termination Type	Moulded-on type K ANSI miniature thermocouple plug
Max. Temperature	260°C
Min. Temperature	-75°C
Tolerance	Tolerance Class 1 for optimum accuracy
Standards Met	ANSI MC96.1 Colour Coded