



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC-Luftsonden-Thermoelement-Sonde mit angespritztem Miniaturstecker - Typ T - 110 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

316 Edelstahl-Luftsonde mit belüftetem Mantel für die allgemeine Messung der Lufttemperatur.

Leistungsbeschreibung

Specifications

Product Code	XE-0014-001
General Description	Vented sheath for general air temperature measurement
Sensor Type	Mineral insulated
Thermocouple Type	T
Probe Diameter	4.0mm
Probe Length	110mm
Termination Type	Miniature moulded flat pin plug, colour coded in accordance IEC specification
Max. Temperature	+250°C (Moulded-on plug +110°C)
Standards Met	IEC
Connector Type	Plug
Connector Size	Miniature (moulded on)
Probe Material	vented sheath, 316 stainless steel