



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

IEC Mineralisoliertes Dual-Thermoelement mit Standard-Duplex-Typ-K-Stecker - 6,0 mm x 250 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

Bitte beachten Sie, dass Sie bezüglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.

Leistungsbeschreibung**Specifications**

Product Code	XE-2014-001
General Description	310 stainless steel sheath. Sheath can be bent/formed to suit many applications and processes.
Sensor Type	Mineral insulated Type 'K' with Duplex Thermocouple Plug
Thermocouple Type	K
Thermocouple Junction	Insulated hot junction
Probe Diameter	6.0mm
Probe Length	250mm
Termination Type	Standard duplex plug termination (200°C)
Max. Temperature	1100°C
Min. Temperature	-40°C
Standards Met	Conforms to IEC 584 specification
Connector Type	Plug
Connector Size	Standard