



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Allzweck-Pt100-Sonde mit PFA-Leitung - Typ PRT - 3,0 mm Sondendurchmesser - 25 mm Sondenlnge - 2 m Kabellnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Allzweck-Pt100-Sonde mit PFA-Leitung - Typ RTD / PRT

Mehrzweck-Pt100-Sonde nach IEC751

Klasse B, 4-Draht-Anschluss

Starrer Mantel aus Edelstahl 316

Teflonisoliertes, geschirmtes Blei

Temperaturbereich des Sensors -50 C bis +250 C

Leistungsbeschreibung**Specifications**

Product Code	XE-3605-001
General Description	Multipurpose Pt100 Probe to IEC751
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Sheath Material	3.0mm diameter 316 Stainless Steel Sheath
Cable Length	2M
Cable Type	Flexible Teflon® insulated, screened
Termination Type	4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	IEC751
Series	General Purpose Pt100