

**Allzweck-Pt100-Sonde mit PFA-Leitung - Typ PRT  
- 4,5 mm Sondendurchmesser - 200 mm  
Sondenlnge - 2 m Kabellnge**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Allzweck-Pt100-Sonde mit PFA-Leitung - Typ RTD / PRT**

Mehrzweck-Pt100-Sonde nach IEC751

Klasse B, 4-Draht-Anschluss

Starrer Mantel aus Edelstahl 316

Teflonisoliertes, geschirmtes Blei

Temperaturbereich des Sensors -50 C bis +250 C

**Leistungsbeschreibung****Specifications**

|                            |                                                                       |
|----------------------------|-----------------------------------------------------------------------|
| <b>Product Code</b>        | XE-3621-001                                                           |
| <b>General Description</b> | Multipurpose Pt100 Probe to IEC751                                    |
| <b>Sensor Type</b>         | Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire                    |
| <b>Sheath Material</b>     | 4.5mm diameter 316 Stainless Steel Sheath                             |
| <b>Cable Length</b>        | 2M                                                                    |
| <b>Cable Type</b>          | Flexible Teflon® insulated, screened                                  |
| <b>Termination Type</b>    | 4-wire tails allows connection to any PT100 instrument 2, 3 or 4 wire |
| <b>Max. Temperature</b>    | 250°C                                                                 |
| <b>Min. Temperature</b>    | -50°C                                                                 |
| <b>Standards Met</b>       | IEC751                                                                |