


Handnadelsonde T-förmig - Asphalt- oder Lebensmittelverarbeitung 1M / 1.4M - Typ K



[Click here to view product](#)

Product Code	XE-3798
General Description	probe response time: less than 17 seconds
Sensor Type	T-shaped handle reduced tip penetration probe
Thermocouple Type	K
Probe Diameter	9.5mm
Probe Length	1000 mm 1400 mm
Max. Temperature	250°C
Min. Temperature	-75°C
Standards Met	BS EN 60584-1:2013