



Tragbare Oberflächensonde - Typ K



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mantel aus rostfreiem Stahl mit geformtem Griff und 2 Meter Spiralverlängerungskabel, das in einem Miniatur-Thermoelementstecker endet.

Specifications

Product Code	XE-3840-001
General Description	Handheld Surface Probe
Thermocouple Type	K
Calibration	IEC
Probe Diameter	15mm
Probe Length	170mm
Sheath Material	Stainless steel sheath
Cable Length	2m
Cable Termination	Fitted Plug
Max. Temperature	400°C
Series	Handheld Probe