

## **IEC Krokodilklemmen-Thermoelement mit berflochtenem Glasfaser-Edelstahlkabel - Typ J - 3M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                                                               |
|------------------------------|-------------------------------------------------------------------------------|
| <b>Product Code</b>          | XE-5651-001                                                                   |
| <b>General Description</b>   | Crocodile Clip Thermocouple with Glassfibre Stainless Steel Overbraided Cable |
| <b>Thermocouple Type</b>     | J                                                                             |
| <b>Thermocouple Junction</b> | Grounded measuring junction                                                   |
| <b>Cable Length</b>          | 3M                                                                            |
| <b>Cable Type</b>            | Stainless steel over-braided Glassfibre                                       |
| <b>Cable Termination</b>     | Miniature Plug                                                                |
| <b>Max. Temperature</b>      | 250°C                                                                         |
| <b>Min. Temperature</b>      | -50°C                                                                         |
| <b>Standards Met</b>         | IEC                                                                           |
| <b>Series</b>                | Crocodile Clip Thermocouple                                                   |