



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Pt1000 Klasse B Sensoren mit teflonisolierter Leitung - Typ PRT - 4,0 mm Durchmesser - 125 mm Länge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung

Specifications

Product Code	XE-3684-001
General Description	Insulated element in a 4.0mm diameter 316 stainless steel probe
Sensor Type	Pt1000 (1000 Ohms @ 0°C) to IEC 751, Class B
Element	Pt1000
Probe Diameter	4.0mm
Probe Length	125mm
Sheath Material	316 stainless steel probe
Cable Length	1M
Cable Type	Teflon® insulated twisted leads
Cable Termination	Bare tails
Max. Temperature	200°C
Min. Temperature	-50°C
Number of Wires	2
Series	Pt1000 Sensors