



Miniatur-Thermoelement-Stecker IM-K-M Typ K IEC



Miniature in-line Plug
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Stecker IM-K-M Typ K IEC

Mit flachen Pins passend zum Thermoelementtyp und ideal für den allgemeinen Einsatz sind alle Kontakte polarisiert, um eine korrekte Verbindung zu gewährleisten. Diese Steckverbinder können Thermoelementkabel mit bis zu 4 mm Durchmesser aufnehmen.

Spezifikationen

Specifications

| | |
|----------------------------|---|
| Product Code | XE-1100-001 |
| General Description | Ideal for general purpose use. Contacts are polarised to ensure correct connection. |
| Thermocouple Type | K IEC |
| Connector Type | Plug |
| Connector Size | Miniature |
| Colour | Green |
| Max. Temperature | 220°C |
| Pins | Flat Pins |
| +Positive Contact | Nickel Chromium |
| -Negative Contact | Nickel Alloy |
| Standards Met | BS EN 50212:1997. RoHS |
| Series | Miniature Connector |