


## Miniatur-Thermoelement-Steckverbinder In-Line-Buchse IM-K-FS Typ K IEC



**Miniature in-line Socket**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur-Thermoelement-Steckverbinder In-Line-Buchse IM-K-FS Typ K IEC**

Ideal for general use, all contacts are polarised, to ensure correct connection.  
These connectors can accept thermocouple cables with a diameter of up to 4 mm.

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1102-001
<b>General Description</b>	Ideal for general purpose use. Contacts are polarised to ensure correct connection.
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector