


## Blaues Klappfhler-Thermometer



[Click here to view product](#)

<b>Product Code</b>	XE-3952-001
<b>Model</b>	TP305
<b>General Description</b>	Household, folding thermometer, digital folding BBQ grill thermometer.
<b>Type</b>	Celsius and Fahrenheit
<b>Accuracy</b>	+/-1°C at -20°C to 200°C. +/-2°C at other ranges
<b>Resolution</b>	0.1 C/F
<b>Response</b>	Sampling time: 1 second / Response time: 4 seconds
<b>Battery</b>	AAA battery powered (included)
<b>Dimensions</b>	270mm x 35mm x 22mm. Probe length: 120mm / Probe diameter: 3mm
<b>Functions</b>	Fold / unfold the probe to switch on / off

<b>Function #1</b>	10 minutes auto-off
<b>Function #2</b>	Super fast response time
<b>Probe Diameter</b>	3mm
<b>Probe Length</b>	120mm
<b>Probe Material</b>	ABS and stainless steel
<b>Sensor Type</b>	Temperature Sensor
<b>Sheath Material</b>	ABS plastic case, LCD digital Display
<b>Ranges Â°C</b>	-50°C to 300°C