



Miniatur Quick Wire Connector Thermoelement-Stecker FMTC-CU-MQ Typ Cu



Miniature in-line Plug - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur-Thermoelement-Steckverbinder Schnelldrahtstecker FMTC-CU-MQ Typ CU

Miniaturstecker mit flachen Steckerstiften. Einfache Jab-In-Verbindung fr eine schnelle Kndigung. Einfach den Draht eindrcken und die Schraube festziehen.

Technische Daten**Specifications**

Product Code	XE-1381-001
General Description	Easy jab-in connection for rapid termination.
Thermocouple Type	IEC, ANSI, JIS Type CU
Connector Type	Quick Wire Plug
Connector Size	Miniature
Colour	White
Max. Temperature	220°C
Pins	Flat Pins
+Positive Contact	Copper
-Negative Contact	Copper
Standards Met	BS EN 50212:1997. RoHS