


Miniatür Quick Wire Connector Thermoelement-Buchse IM-K-FQ Typ K IEC



Miniature in-line Socket - Quick Wire
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniatur Quick Wire Connector Thermoelement-Buchse IM-K-FQ Typ K IEC

Einfache Einstichverbindung, einfach den Draht eindrücken und die Schraube festziehen. Die Kontakte sind polarisiert, um eine falsche Terminierung zu verhindern.

Technische Daten**Specifications**

Product Code	XE-1112-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	K IEC
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS