


## Verstellbares Bajonett, berflochtenes Glasfaserkabel aus Edelstahl - Typ PRT - 1M Lnge



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Verstellbares Bajonett mit 200 mm Feder mit verstellbarer Befestigungskappe (11,5 mm Innendurchmesser Einzelschlitz).

### Leistungsbeschreibung

## Specifications

<b>Product Code</b>	XE-5654-001
<b>General Description</b>	PRT with adjustable insertion depth
<b>Sensor Type</b>	Pt100 4-wire
<b>Cable Length</b>	1M
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel Overbraid
<b>Core / Strands</b>	4 core 7/0.2mm
<b>Cable Termination</b>	Bare tails
<b>Max. Temperature</b>	350°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	IEC-60751 Class B
<b>Bayonet cap</b>	11.5mm I.D. single slot
<b>Spring</b>	6.0mm diameter x 200mm long