

**IEC Typ K 9kg Pull Power (Hufeisen)
Magnet-Thermoelement, PFA-isoliertes Kabel mit
Edelstahl-bergeflecht mit Abschluss im
Miniaturstecker - 3,0 m**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-5641-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Sensor Type	Type K thermocouple to IEC 584
Thermocouple Type	K
Thermocouple Junction	Grounded at tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Length	3.0M
Cable Termination	PFA Teflon® insulated with stainless steel over-braid, miniature plug termination
Termination Type	Miniature Plug
Max. Temperature	+ 250°C
Min. Temperature	-50°C
Standards Met	IEC