


# Miniatür Quick Wire Connector Thermoelement-Buchse FMTC-CU-FQ Typ Cu



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Quick Wire Connector Thermoelement-Buchse FMTC-CU-FQ Typ CU**

Einfache Einstichverbindung, einfach den Draht eindrücken und die Schraube festziehen. Die Kontakte sind polarisiert, um eine falsche Terminierung zu verhindern.

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1382-001
<b>General Description</b>	Easy jab-in connection, just push in wire & tighten screw.
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Quick Wire Socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS