

## **IEC Typ K 3kg Zugmagnet-Thermoelement, PFA-isoliertes Kabel mit bergeflochtenem Edelstahlgeflecht mit Miniaturstecker - 1M**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                                         |
|------------------------------|---------------------------------------------------------|
| <b>Product Code</b>          | XE-5616-001                                             |
| <b>General Description</b>   | Rapid, reliable temperature sensing on ferrous surfaces |
| <b>Thermocouple Type</b>     | K                                                       |
| <b>Thermocouple Junction</b> | Grounded measuring junction                             |
| <b>Magnet</b>                | 1" diameter (3kg pull)                                  |
| <b>Cable Length</b>          | 1M                                                      |
| <b>Cable Type</b>            | Stainless steel over-braided PFA Teflon®                |
| <b>Cable Termination</b>     | Miniature Plug                                          |
| <b>Max. Temperature</b>      | +250°C (200°C Plug)                                     |
| <b>Standards Met</b>         | IEC                                                     |