




LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**IEC-Unterlegscheiben-Thermoelement,
glasfaserisoliertes Kabel mit
Edelstahl-bergeflecht Typ K mit blanken
Endverschlüssen - 4,0 mm
Unterlegscheibenlochgröße**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Leistungsbeschreibung**Specifications**

Product Code	XE-3428-001
General Description	Thermocouple housed in a washer eyelet
Thermocouple Type	K
Thermocouple Junction	Grounded tip for quick response
Cable Length	2m
Cable Type	Glassfibre insulated, stainless steel overbraided
Core / Strands	7/0.2mm
Termination Type	Bare Tails
Max. Temperature	+350°C
Standards Met	IEC
Washer	3.5mm, 4.0mm, 5.0mm, 6.0mm, 8.0mm & 10.0mm clearance hole available
Series	Washer Thermocouple