




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

**IEC-Unterlegscheiben-Thermoelement,  
glasfaserisoliertes Kabel mit  
Edelstahl-bergeflecht Typ K mit blanken  
Endverschlüssen - 6,0 mm  
Unterlegscheibenlochgröße**



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Leistungsbeschreibung****Specifications**

|                              |                                                                     |
|------------------------------|---------------------------------------------------------------------|
| <b>Product Code</b>          | XE-3416-001                                                         |
| <b>General Description</b>   | Thermocouple housed in a washer eyelet                              |
| <b>Thermocouple Type</b>     | K                                                                   |
| <b>Thermocouple Junction</b> | Grounded tip for quick response                                     |
| <b>Cable Length</b>          | 2m                                                                  |
| <b>Cable Type</b>            | Glassfibre insulated, stainless steel overbraided                   |
| <b>Core / Strands</b>        | 7/0.2mm                                                             |
| <b>Termination Type</b>      | Bare Tails                                                          |
| <b>Max. Temperature</b>      | +350°C                                                              |
| <b>Standards Met</b>         | IEC                                                                 |
| <b>Washer</b>                | 3.5mm, 4.0mm, 5.0mm, 6.0mm, 8.0mm & 10.0mm clearance hole available |