


Miniaturrechteckiger Thermoelement-Steckverbinder Faszienbuchse IM-T-FF Typ T IEC



Miniature Fascia Socket (with nylon mounting clip)
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniaturrechteckiger Thermoelement-Steckverbinder Faszienbuchse IM-T-FF Typ T IEC

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

Technische Daten
Specifications

Product Code	XE-1197-001
General Description	All fascia mounting sockets are polarised by means of unequally sized pins
Thermocouple Type	T IEC
Connector Type	Fascia mounting socket
Connector Size	Miniature
Colour	Brown
Dimensions)
Max. Temperature	120°C
+Positive Contact	Copper
-Negative Contact	Copper Nickel
Standards Met	BS EN 50212:1997. RoHS
Series	Miniature Connector