



# Miniaturrechteckiger Thermoelement-Steckverbinder Faszienfassung FMTC-CU-FF Typ Cu



**Miniature Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniaturrechteckiges Thermoelement-Steckverbinder-Faszien-Buchse FMTC-CU-FF Typ CU**

Nylonclip für schnelle, einfache Befestigung am Panel. Alle Fassadenmontage-Fassungen sind durch ungleich große Pins polarisiert

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1377-001
<b>General Description</b>	All fascia mounting sockets are polarised by means of unequally sized pins
<b>Thermocouple Type</b>	IEC, ANSI, JIS Type CU
<b>Connector Type</b>	Fascia mounting socket
<b>Connector Size</b>	Miniature
<b>Colour</b>	White
<b>Max. Temperature</b>	120°C
<b>+Positive Contact</b>	Copper
<b>-Negative Contact</b>	Copper
<b>Standards Met</b>	BS EN 50212:1997. RoHS