


# Miniatur Quick Wire Thermoelement Steckverbinder-Montage mit Edelstahlhalterung IM-K-SSPFQ Typ K IEC



**Miniature Socket -**

**Quick Wire with stainless steel mounting bracket**

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniatur Quick Wire Thermoelement Steckverbinder-Montage mit Edelstahlhalterung IM-K-SSPFQ Typ K IEC**

Miniatur-Fassungen mit Edelstahlhalterungen ermöglichen eine schnelle, zuverlässige und robuste Montage in das Panel

**Technische Daten****Specifications**

<b>Product Code</b>	XE-1106-001
<b>General Description</b>	Contacts are polarized to ensure correct connection
<b>Thermocouple Type</b>	K IEC
<b>Connector Type</b>	Quick Wire Socket with SS Bracket
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Max. Temperature</b>	220°C
<b>+Positive Contact</b>	Nickel Chromium
<b>-Negative Contact</b>	Nickel Alloy
<b>Standards Met</b>	BS EN 50212:1997. RoHS
<b>Series</b>	Miniature Connector