



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Mineralisiertes Thermoelement mit 0,25 mm Durchmesser und schnellem Ansprechverhalten mit IEC-Miniaturstecker - Typ K - geerdete Verbindungsstelle - 100 mm



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Mssen Sie Ihren Temperatursensor kalibrieren?

Labfacility ist auf ISO17025 (UKAS)-Kalibrierstandards rckfhrbar und bietet die Wahl zwischen einer 3- oder 5-Punkt-rckfhrbaren Kalibrierung mit unserer hauseigenen Kalibriereinrichtung.

Bitte beachten Sie, dass Sie bezüglich der Lieferinformationen kontaktiert werden, wenn Sie Ihrer Bestellung eine Kalibrierungsoption hinzufügen.

Leistungsbeschreibung**Specifications**

Product Code	XE-0806-001
General Description	Highly flexible, sheath can be bent/formed to suit many applications and processes
Sensor Type	Type K IEC MB STYLE probes 0.25mm
Thermocouple Junction	Single element, grounded junction
Probe Length	100mm
Sheath Material	Stainless steel 316 sheath
Termination Type	Fitted Miniature Plug Type K IEC
Max. Temperature	260°C
Min. Temperature	-40°C
Important Information	A 1%+ Length Tolerance applies