

Lascar EL-P-TPX+ kalibrierbare Thermistor-Sonde



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Hochprzise kalibrierbare, 3 Meter lange Thermistorsonde mit Edelstahlspitze fr den Einsatz mit EL-WiFi-TPX+ ([XE-4076-001](#))

Technische Daten**Specifications**

Product Code	XE-4081-001
Model	EL-P-TPX+
Type	Calibratable thermistor probe
Accuracy	+/-0.2°C (-15 to +80°C) / +/-0.6°C (-40 to 125°C)
Cable Type	3m (118.1")
Measuring range	-40 to +125°C / -40 to + 257°F